

5.1 Materials provided

1. Lipase substrate and buffer reagents.
2. Lipase activator solution.
3. Lipase standard.

<u>Mean</u>	<u>S.D.</u>	<u>C.V.%</u>
46.7	1.70	3.64
254.0	1.70	1.47
516.5	4.65	0.90

5.2 Materials Required but not Provided

1. Accurate pipettes.
2. Test tubes.
3. Timer.
4. Heating block or water bath (37°C)
5. Spectrophotometer with a temperature controlled cuvette.
6. Controls

6.0 Test Procedure

6.1 Procedure (Automated)

See appropriate instrument application instructions.

6.2 Procedure (Manual)

1. Label test tubes Blank, Standard, Control, Patient, etc.
2. Pipette 300µl of reconstituted Lipase Substrate reagent to all tubes.
3. Pipette 5µl of distilled water to the blank tube and 5µl of the appropriate sample to the tubes labeled Standard, Control, etc.
4. Mix each tube well and incubate for 3-5 minutes at 37°C.
5. After the pre-incubation, add 100µl of Lipase activator to the blank tube. Mix well and incubate for 3 minutes at 37°C. Then measure the rate of increase in absorbance per minute at 550nm (540-560nm).
6. Repeat step 5 for all tubes.
7. See Calculations to obtain results.

6.3 Procedure Notes

The above volumes may be multiplied by an appropriate factor if a larger total volume of reaction mixture is necessary for reading.

6.4 Calibration

Use the lipase standard provided in the kit.

6.5 Quality Control

The integrity of the reaction should be monitored by use of normal and abnormal control sera with known lipase activity.

6.6 Limitations

Samples with Lipase activity exceeding 600 U/L should be diluted with an appropriate amount of saline, re-assayed, and the final result multiplied by the appropriate dilution factor.

6.7 Calculation

$$\frac{\Delta A \text{ Sample} - \Delta A \text{ Blank}}{\Delta A \text{ Standard} - \Delta A \text{ Blank}} \times \text{Conc of Std. (U/L)} = \text{Lipase activity (U/L)}$$

6.8 Expected Values

0-62 U/L

It is strongly recommended that each laboratory establish its own normal range.

7.0 Performance

1. Linearity: 600U/L
2. Comparison: A study performed comparing the Lipase (Colorimetric) methods to a turbidimetric Lipase procedure yielded a correlation coefficient of 0.956 with a regression equation of $y = 0.48x + 9.1$.
3. Precision: The lipase activity of three samples was measured ten times each with the following results:

8.0 References

1. Imamura, S. and Misaki, H., Selected Topics in Clinical Enzymology, 2:73, 1984.
2. Imamura S. et al, Collection of summaries of lectures in the 126th general meeting of Kinki branch, analytical section, Japan Society of Clinical Chemistry, p 11-31, 1986.
3. Hayashi, C., et al, Clinical Examination, Instrument and Reagent, 2:25, 1986.
4. Kitaura, S., et al, Collection of summaries of lectures in the 6th general meeting of The Japanese Biochemical Society, 848, 1988.
5. Imamura, S., et al, clin., chem., Abstract Issue in the 41st National Meeting, 1120, 1989.
6. Young, D.S., et al, Clin. Chem. 21:1D, 1975.